


EW-Type

Impression Parts

Impression Parts

Common to BW-type




Set name	W-ME Coping-S impression	W-M3 Coping-S impression	N-M3 Coping-L impression
Impression parts	W-ME Coping-S	W-M3 Coping-S	W-M3 Coping-L
Pin	Incl.	Incl.	Incl.
Prosthetic parts	ReplicaBW-HEX		


Open-tray technique
Fixture Level Impression

Closed-tray technique
Fixture Level Impression

Closed-tray technique
Fixture Level Impression



EW60**H
ReplicaBW60-HEX




Set name	BW-Screw impression
Impression parts	BW-Coping
Pin	L-pin
Prosthetic parts	ReplicaBW-HEX

Open-tray technique
Fixture Level Impression

※For BW-Screw impression only EW50**H is used.

※ Each single parts are available.


Full Set (Indicated for Single Tooth Restorations) For EW50**H only



Set name	AB43-G35(fit)	AB43-G50(fit)	AB43-G70(fit)
Abutments	AB43-G35	AB43-G50	AB43-G70
Prosthetic parts	ReplicaBW-35G	ReplicaBW-50G	ReplicaBW-70G
Impression parts	HAT43G (Common parts)		
Temporary Crown parts	Pattern43G (Common parts)		
Wax-up parts	Pattern43G (Common parts)		

Single Tooth Restorations Closed-tray technique Abutment Level Impression
Single parts are available.

Full Set (Indicated for Multiple-unit Restorations) For EW50**H only



Set name	AB43-J35(fit)	AB43-J50(fit)	AB43-J70(fit)
Abutments	AB43-J35	AB43-J50	AB43-J70
Prosthetic parts	ReplicaBW-35J	ReplicaBW-50J	ReplicaBW-70J
Impression parts	HAT43J (Common parts)		
Temporary Crown parts	Pattern43J (Common parts)		
Wax-up parts	Pattern43J (Common parts)		

Multiple-unit Restoration Closed-tray technique Abutment Level Impression
Single parts are available.

<Notice>

For EW60**H, use the ReplicaBW60-HEX.

For EW70**H, the Replica is not available now. Take the impression by an direct impression.