


























Tools & Drills

Drills

Model name		Full length/Diameter
Round Bur Round Bur	 Used to mark the location of drilling.	Full length:34mm / ϕ 2.0mm
Lindemann Bone Cutter Lindemann Bone Cutter (164RF/018)	 Used to drill cortical bone.	Full length:34mm / ϕ 1.7mm
Lancet Drill Lancet Drill	 Used to drill cortical bone.	Full length:34mm / ϕ 1.67mm
Guide Drills Spiral Drill 16S		Full length:32mm / ϕ 1.6mm
Spiral Drill 20S	 Determine the direction and depth of the implant slot according to the length of the endosteal part of the implant.	Full length:32mm / ϕ 2.0mm
Drills Spiral Drill 25S		Full length:32mm / ϕ 2.5mm
Spiral Drill 28S		Full length:32mm / ϕ 2.8mm
Spiral Drill 32S		Full length:32mm / ϕ 3.2mm
Spiral Drill 35S		Full length:32mm / ϕ 3.5mm
Spiral Drill 38S		Full length:32mm / ϕ 3.8mm
Spiral Drill 41S		Full length:32mm / ϕ 4.1mm
Spiral Drill 45S		Full length:32mm / ϕ 4.5mm
Spiral Drill 47.5S		Full length:32mm / ϕ 4.75mm

Model name		Full length/Diameter
Trephine Burs Hole Saw 5.35		Full length:35mm ϕ 5.35mm (External diameter) ϕ 4.35mm (Internal diameter)
Hole Saw 6.35		Full length:35mm ϕ 6.35mm (External diameter) ϕ 5.35mm (Internal diameter)
End Mills End Mill 5.55		Full length:34mm ϕ 5.55mm
End Mill 6.55		Full length:34mm ϕ 6.55mm
Countersinks Countersink 40		Full length:34.5mm ϕ 4.0mm (Maximum diameter)
Countersink 50		Full length:34.5mm ϕ 5.0mm (Maximum diameter)
Over Drills Over Drill 10J		Full length:26mm ϕ 1.0mm (Cuff)
Over Drill 20J		Full length:26mm ϕ 2.0mm (Cuff)
Over Drill 30J		Full length:26mm ϕ 3.0mm (Cuff)
Circular Knives Circular knife 40		Full length:30mm ϕ 4.0mm
Circular knife 48		Full length:30mm ϕ 4.8mm
Drill Stoppers Drill Stopper		Full length:8.5mm