

Multi-Piece (MP) Abutments for Screw-Retained Restorations

NEW




Multi-Piece (MP) Abutment for Screw-Retained Restorations

【Normal-Connection type】


Compliant implants	C-type, E-type, St-type, Me-type, B-type
--------------------	--

Taper Connection type



Cuff Length	Product name
1.0	MP-AB35C10-s
2.0	MP-AB35C20-s
3.0	MP-AB35C30-s
4.0	MP-AB35C40-s
5.0	MP-AB35C50-s

Internal HEX Connection type



Cuff Length	Product name
1.0	MP-AB35C10-HEX
	An abutment screw is included.
2.0	MP-AB35C20-HEX
	An abutment screw is included.
3.0	MP-AB35C30-HEX
	An abutment screw is included.
4.0	MP-AB35C40-HEX
	An abutment screw is included.
5.0	MP-AB35C50-HEX
	An abutment screw is included.

17° Angled-type




Accessory: MP-AG Screw

Cuff Length	Product name
2.5	MP-AG35C25-HEX17
3.5	MP-AG35C35-HEX17
4.5	MP-AG35C45-HEX17
5.5	MP-AG35C55-HEX17

※ A MP-AG Screw is included.

30° Angled-type




Accessory: MP-AG Screw

Cuff Length	Product name
3.5	MP-AG35C35-HEX30
4.5	MP-AG35C45-HEX30

※ A MP-AG Screw is included.

Cylinder for Multi-Piece



Accessory: MP-Cylinder Screw

Size	Product name
Φ3.3	MP-Cylinder-33
Φ4.3	MP-Cylinder-43

※ A MP-Cylinder Screw is included.

※ Up to 6mm from the cylinder top, you can grind.

Caps



Size	Product name
4.7	MP-CAP-47
6.7	MP-CAP-67

MP Pattern



MP Replica



Kinds	Product name
Lab. Pattern	MP Pattern
Lab. Replica	MP Replica
Impression Coping	MP Coping-OPT
Lab. Screw	Lab. MP-Cylinder Screw (2pcs)
	Lab. MP-AG Screw (2 pcs)
Abutment Carry Tool	MP Holder
Tightening Driver	CA-MP Driver
	WR-MP Driver

MP Coping-OPT



Lab. Screw



Lab. MP-Cylinder Screw




Lab. MP-AG Screw

MP Holder



A carry of a MP-AG Abutment

CA-MP Driver



Compliant abutments: Taper Connection & Internal HEX Connection straight-type

WR-MP Driver



Compliant abutments: Taper Connection & Internal HEX Connection straight-type

[Abutment Tightening Torque Value]

- Taper Connection type
- Internal HEX Connection type



30Ncm



15Ncm



Be sure to use a Cylinder for fabrication of super-structures.

- Angled-type 17°, 30°



30Ncm



15Ncm



10Ncm



※ Use a HEX Driver

[Cap Tightening Torque Value]

Be sure to tighten a MP-Cylinder Screw with less than 15Ncm.

MYTIS Arrow Implant

BrainBase Corporation

Access Oimachi Bldg 6F, 49-15, Oi 1-chome, Shinagawa-ku Tokyo, Japan 140-0014

TEL : +81-3-3778-0745 FAX : +81-3-3778-4910

E-mail : mail@brain-base.com URL <http://www.brain-base.com>

